


FIBION MAINTENANCE GUIDE

ENSURE A LONG LIFE FOR YOUR FIBION


1.



1. Keep it clean and dry

Regularly clean your Fibion Device, especially after working out and sweating. To clean your Fibion Device, do not use household cleansers. Instead, use a soft towel. For tough spots or stains, scrub with a wet, soft-bristled toothbrush. The Fibion Device does not like water. If your Fibion Device gets wet—like after sweating—dry it thoroughly before putting it back in your pocket or thigh strap.

2.



2. Give your Leg a Rest when Using the Thigh strap

When you are using the thigh strap for a longer period of time, ensure that the band isn't too tight. Prolonged rubbing and pressure may irritate the skin, so give your leg a break by removing the thigh strap during the night, unless you are instructed otherwise by a medical trained professional. If you have eczema, allergies or asthma, you may be more likely to experience a skin irritation or allergy from the thigh strap. Whether you have the conditions above or not, if you start to experience redness or skin irritation on your leg, remove your Fibion Device. If symptoms persist longer than two to three days of not using your Fibion Device, contact a dermatologist.

3.



3. Wash the Thigh Strap and Pouches regularly

Hand wash the strap and pouches in lukewarm water. Use only mild detergent. Rinse thoroughly and drip dry only.